

DyNAtrix®

Troubleshooting Guide

Ready-to-Use Kit (#DM-RUK-1) and Minimal Discovery Kit (#DM-MIDK-1) | Version 1.1

1.	Before cell embedding	2
2.	During cell embedding	2
3.	After cell embedding & in culture	2
4.	Cell release from DyNAtrix.....	3
5.	Immunofluorescent staining in DyNAtrix.....	4
6.	Information to provide when contacting technical support.....	4

1. Before cell embedding

Issue	Possible causes	Recommendations
Component is yellow or turbid	• Possible microbial contamination. • pH change caused by improper storage.	• Do not use the component; contact support with the lot number and photographs.
Component is unusually viscous, partially gelled, or difficult to pipette	• Cross-contamination between components initiated partial crosslinking.	• Use a new sterile tip for each stock vial. • Contact support if the stock vials are partially gelled.

2. During cell embedding

Issue	Possible causes	Recommendations
Matrix remains liquid or does not gel after incubation for 60 min at 37 °C	• The final crosslinker was omitted, incorrect, or mismatched. • Serum or conditioned medium was carried into the matrix.	• Verify component volume and the matched crosslinker pair for the Minimal Discovery Kit. • Wash cells twice in serum-free medium before embedding to minimize serum carryover.
Gelation happens before dispensing	• Components or the mixture warmed before dispensing. • Final crosslinker was added to the mixture and sat too long before dispensing.	• Keep components and the mixture cold at 4 °C or on ice during preparation. • Add the final crosslinker immediately before dispensing.
Bubbles are introduced during mixing or dispensing	• Mixing was too rapid, or the full volume was expelled forcefully with air mixed into the mixture.	• Keep the tip below the liquid surface and pipette about 80% of the total volume. • Quickly spin the tube on a mini benchtop centrifuge to remove bubbles. • Remove a surface bubble carefully after dispensing with a sterile tip without disturbing the gel.
Cell-matrix mixture is too viscous or difficult to pipette	• B+ is slightly viscous by design. • Partial gelation has started.	• Use a positive displacement pipette to mix and dispense for accurate volume and low material retention in tips.
Cells are unevenly distributed within gels	• Incomplete mixing. • The cell stock settled or contained clumps before the aliquot was taken.	• Mix by pipetting up and down 20 µL for 25 times for homogenous mixing; avoid bubbles. • Gently resuspend the cell stock before taking the aliquot.
Spheroids or organoids look fragmented after mixing or dispensing	• The pipette-tip opening was too narrow, causing excessive shear. • Spheroids were too large.	• Use a wide-bore tip for mixing and dispensing. • Prepare spheroids no larger than 100 µm in diameter before embedding to reduce shear stress.
Mixture looks inhomogeneous or phase-separated	• Incomplete mixing, or an incorrect order of addition.	• Prepare a fresh sample using the exact volumes and order; mix thoroughly by pipetting up and down for ≥25 times.

3. After cell embedding & in culture

Issue	Possible causes	Recommendations
Gel dome detaches or is lost during medium handling	• Medium was added or removed directly onto the gel.	• Add and aspirate medium along the vessel wall. Use a microwell insert for small gels or drop a larger dome volume.

Issue	Possible causes	Recommendations
Dome breaks or collapses	Incomplete mixing or gelation.	- Prepare a fresh sample using the exact component volumes and order; mix thoroughly by pipetting up and down for ≥ 25 times. - Incubate the mixture at 37 °C for 60 min to allow complete gelation.
Gel degrades unintentionally, and cells sediment to the bottom	- Serum, KnockOut Serum Replacement, conditioned medium, or another nuclease-containing supplement was used. - Serum or conditioned medium was carried in with the cells. - Extensive cell death released nucleases.	- Use fully defined, serum-free medium throughout the culture unless a formulation has been specifically validated. - Wash cells twice in serum-free medium before embedding. - Contact support with the medium composition if degradation persists.
Spheroids or organoids sediment to the bottom during gelation	Large or dense aggregates weight heavier and sediment before gelation is complete.	Flip the plate upside down for the first 10 min, then return upright for the rest 50 min of incubation at 37 °C.
Little or no cell proliferation or organoid expansion	- The defined medium lacks cell-type-specific growth factors or nutrients. - The RGD ligand in the matrix is insufficient for the model.	- Run a cell-only control on 2D or in suspension to confirm viability in the serum-free medium. - Supplement B+ with 0.2-0.5 mg/mL (recombinant) laminin before gelation. - Contact support to discuss alternative adhesion ligands.
Unhealthy morphology or cell death	- Excessive pipetting shear occurred during mixing or dispensing. - The defined medium is not optimal for the cell model. - The RGD ligand in the matrix is insufficient for the model.	- Use wide-bore tips for cell embedding. - Contact support for serum-free medium suggestions. - Contact support to discuss alternative adhesion ligands.

4. Cell release from DyNAtrix

Issue	Possible causes	Recommendations
Gel degradation is slow or incomplete after DNase I treatment	DNase I concentration or activity was insufficient.	- Add more DNase I to the medium or extend incubation time. Use fresh, active DNase I for optimal performance. - Gently pipette the gels to mechanically break them and facilitate degradation.
Released cells adhere to the vessel before collection	The gel was digested in a tissue-culture-treated vessel that promotes cell adhesion.	- Transfer the gel to a low-binding tube containing 0.2-0.5 mL medium before adding DNase I. - Pre-wet the tube and tip with 1% BSA solution to cover the surface with anti-sticking coating.
Residual gel or DNase I interferes with downstream analysis	Incomplete gel degradation or insufficient washing.	Wash released cells by centrifuging and resuspending the pellet in fresh, assay-compatible buffer or serum-free medium. Verify that visible gel fragments are gone.

5. Immunofluorescent staining in DyNAtrix

Issue	Possible causes	Recommendations
DAPI or Hoechst nuclear signal is weak in cells or the matrix background is strong	DyNAtrix contains DNA-based crosslinks that can be stained by DNA-binding nuclear dyes.	- Increase dye concentration for ≥ 5 folds. - Release the cells before staining.
Antibody signal is weak or visible only near the gel surface	Antibody penetration is slow in thick 3D samples.	- Extend antibody incubation time and gently rotate the plate while keeping the gel fully immersed. - Use Fab fragments or nanobodies instead of full-length antibodies for faster penetration. - Use a smaller or thinner gel format. - Release the cells before staining.
High nonspecific binding of antibody in matrix background	- Washing was insufficient. - Unbound antibody did not diffuse out of the gel.	- Incubate in the washing solution (1%BSA in PBS) overnight after primary antibody staining to remove unbound antibody. Gently rotate the plate while keeping the gel fully immersed. - Repeat extended washing by incubating in the washing solution overnight after secondary-antibody staining. - Use Fab fragments or nanobodies instead of full-length antibodies for faster diffusion.
Gel detaches or floats during imaging	The culture vessel does not retain the sample during repeated washing.	Remove most buffer and cover the floating gel with a coverslip to stabilize its position. Do not allow the gel to dry.
Deep structures are blurred or difficult to image	- Light scattering and refractive-index mismatch in a thick sample. - The culture vessel is not glass-bottom.	- Immerse the gel in mounting media with a refractive index of 1.42-1.45 to minimize light scattering. - Perform optical clearing in DyNAtrix after fixation, such as delipidation in SDS buffer, to increase transparency of large multicellular structures. - Use an imaging-grade vessel with a #1.5 glass bottom or equivalent.

6. Information to provide when contacting technical support

- ☐ Kit name and lot number.
- ☐ Cell type and medium type.
- ☐ Culture vessel type.
- ☐ Culture duration.
- ☐ Dyes, antibodies, and staining methods, if relevant.
- ☐ If available, provide photograph or microscopy images to illustrate the issue.