

Materials included in the package

Component	Amount (µL)
Backbone (B+)	1200
X-linker (X)	400

General remarks

This guide describes standard cell culture in dome or microwell format*. The protocol can be adapted for different cell densities or culture formats.

!!! DyNAtrix is optimized for fully defined media. **Do not add serum or conditioned media**, as they may degrade the matrix.


* e.g. ibidi Cat. # 80409

Preparation

1. Prepare cell suspension in a defined, serum-free medium.


Recommended cell-stock density: $1-4 \times 10^6$ cells/mL

2. Pipette 15 µL of B+ into a tube. Keep it on ice.



Cell embedding at 4 °C / on ice

3. Add 5 µL cell suspension to the 15 µL B+.
Mix *thoroughly* by pipetting up and down.
4. Add 5 µL X to the mixture from step 1.
Mix *thoroughly* by pipetting up and down.
5. Immediately dispense 5 µL per microwell or
≥20 µL per dome.



Initiate culture conditions at 37 °C

1. Incubate at 37 °C, 5% CO₂ for 60 min.
2. Gently cover the microwell insert or dome with serum-free medium.
3. Continue incubation at 37 °C, 5% CO₂.



For detailed protocol, visit
www.dynamicmatrices.eu